

Also available from BioPoint:-

Cryptosporidium and Giardia Testing

ColorSeed™

EasyStain™

www.biopoint.com.au

© Copyright 2019 BioPoint Pty Ltd Sydney Australia

Safety Precautions and Limitations

After use, dispose of packaging in accordance with appropriate biohazard disposal practices. Do not use the product if the packaging is damaged. Do not use the product after the expiry date. The product should be used according to the procedure described in this instruction for use document. Any modifications may affect the results.

Warranty

BioPoint warrants to the original purchaser only that the products will conform to the quantity and contents stated on the product labels and quality certificate for the duration of the stated shelf life.

Limits of Liability

BioPoint's sole obligation and the purchaser's exclusive remedy under the above warranty is limited either to replacement, at BioPoint's expense, any products which are defective in manufacture (provided that the purchaser must return the products, transportation prepaid, to BioPoint) or, at BioPoint's option, refund the purchase price. All claims must be notified to BioPoint within 14 days of using the product. The above warranty does not apply to any products which have been altered outside BioPoint, nor to any products which have been subjected to misuse or mishandling, or use, storage or handling other than in accordance with instructions supplied by BioPoint. ALL OTHER WARRANTIES, EXPRESS, IMPLIED OR STATUTORY, ARE HEREBY SPECIFICALLY EXCLUDED, INCLUDING BUT

NOT LIMITED TO WARRANTIES OR MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. BioPoint's maximum liability is limited in all events to the price of the products sold by BioPoint. IN NO EVENT SHALL BioPoint BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES OR LOSS OF PROFITS. The limitations set out above apply to the fullest extent permitted by law. If any provision of these terms and conditions is held to be invalid, illegal or unenforceable, the enforceability of the remaining provisions will be unaffected.

Use and Indemnity

The purchaser assumes all risk and responsibility in connection with the use, storage, handling and disposal of the products. The purchaser must ensure that all use, storage, handling and disposal of the products is in compliance with all applicable laws, regulations and guidelines and in accordance with any instructions supplied by BioPoint, and must take all appropriate safety precautions. The purchaser indemnifies BioPoint and its directors, officers and employees and agents from and against all claims, actions, liabilities, costs and expenses arising directly or indirectly from the use, storage, handling and disposal of the products, except to the extent arising directly from the negligence of BioPoint.

© Copyright 2019 BioPoint Pty Ltd Sydney Australia exercising jurisdiction there, EasySeed™, ColorSeed™, EasyStain™.

EasySeed™
Cryptosporidium Only

Instructions For Use

Issue: 1

Date: October 2019



Instructions for Use

Seeding the sample (use one tube of EasySeed™)

1. Remove and keep the tube cap.
2. Add 2mL of 0.05% (v/v) Tween 20 to the tube.
3. Replace cap and vortex for 20 seconds.
4. Remove and keep cap. Add tube contents to sample.
5. Add 3mL of reagent grade water to the empty tube.
6. Replace cap and vortex for 20 seconds.
7. Remove and keep cap. Add tube contents to sample.
8. Repeat steps 5, 6 and 7 once more.

Sample Analysis

9. Analyze the sample as per the laboratory Standard Operating Procedure.
10. Record the number of *Cryptosporidium* detected.
11. Calculate the *Cryptosporidium* recovery using the following formulae:-

***Cryptosporidium* Recovery (%) =**

$$\frac{\text{Cryptosporidium detected} \times 100}{\text{number of } \textit{Cryptosporidium} \text{ in EasySeed}^{\text{TM}} \text{ as per Certificate of Analysis}}$$

Important

- Do not attempt to dry the contents of EasySeed™ onto a well slide.
- Every Easyseed™ batch is subjected to Quality Control Testing to meet Batch Specifications.

Technical assistance

For technical assistance with the use of EasySeed™ or other *Cryptosporidium* and *Giardia* testing issues please email info@biopoint.com.au

Product Information Sheet

Product Information Sheet can be downloaded from the BioPoint website:
– www.biopoint.com.au

Shipping and Temperature Tolerance

Please refer to the BioPoint Shipping Policy displayed on the BioPoint website
– www.biopoint.com.au