

***Cryptosporidium parvum***  
**10E8**

Product Information Sheet

**Batch number: 740 - 456**

**Production Data**

**Volume:** 1ml  
**Viability:** Viable  
**Catalogue Number:** C10E8-V  
**BTF Production Date:** 29 January 2019

**Stock specifics**

**Strain:** Iowa  
**Source:** Bovine  
**Date purified:** 12 January 2019  
**Date Excreted:** 12 January 2019  
**Purification Method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

**Quality Assurance Specialist**



.....  
**Linh Le**  
4 February 2019



**Water Industry Enquiries:**

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8221 9589  
[www.btfbio.com](http://www.btfbio.com)



Issue: 1 - August 2018