

Batch Number: B 281

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: CSCG100, CS2CG100, CS4CG100, CS50CG100, W-CSCG100, W-CS4CG100
Preparation Date: 22 November 2019
Expiration Date: 21 March 2020
Volume: 1.081 ml +/- 22 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts <small>(Method Ref: CG-011)</small>		
	Mean	SD
Cryptosporidium count:	99	2.1
	Mean	SD
Giardia count:	98	1.5
DAPI staining		
Cryptosporidium % +ve	98	%
Giardia % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 25 October 2019
Purification method: Discontinuous sucrose and cesium chloride centrifugation

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 21 October 2019
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Manager



Lucy Millican
 29 November 2019

Enquiries:

Biopoint Pty Ltd
 Unit 1 35-41 Waterloo Rd
 Macquarie Park NSW 2113
 Tel. 61 (0) 2 8877 9127
www.biopoint.com.au
 info@biopoint.com.au



