

Batch Number: **B 685**

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESCG100, ESCG100-5
Preparation Date: 7 November 2019
Expiration Date: 6 March 2020
Volume: 1.311 ml +/- 11 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts		
<i>(Method Ref: ESCG-001)</i>		
	Mean	St.Dev.
Cryptosporidium count:	99	1.8
	Mean	St.Dev.
Giardia count:	99	1.9
DAPI staining:		
Cryptosporidium % +ve	100	%
Giardia % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 25 September 2019
Purification method: Discontinuous sucrose and cesium chloride centrifugation

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 21 October 2019
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Manager

.....

Lucy Millican

20 November 2019

Water Industry Enquiries:

Biopoint Pty Ltd
 Unit 1 35-41 Waterloo Rd
 Macquarie Park NSW 2113
 Tel. 61 (0) 2 8877 9127
www.biopoint.com.au
 info@biopoint.com.au

