

Batch Number: B 282

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: CSCG100, CS2CG100, CS4CG100, CS50CG100, W-CSCG100, W-CS4CG100
Preparation Date: 12 December 2019
Expiration Date: 10 April 2020
Volume: 1.160 ml +/- 17 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts <small>(Method Ref: CG-011)</small>		
	Mean	SD
Cryptosporidium count:	99	1.6
	Mean	SD
Giardia count:	99	1.7
DAPI staining		
Cryptosporidium % +ve	100	%
Giardia % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 20 November 2019
Purification method: Discontinuous sucrose and cesium chloride centrifugation

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 21 October 2019
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Manager



Lucy Millican
 20 December 2019

Enquiries:

Biopoint Pty Ltd
 Unit 1 35-41 Waterloo Rd
 Macquarie Park NSW 2113
 Tel. 61 (0) 2 8877 9127
www.biopoint.com.au
 info@biopoint.com.au

