

Batch Number: **B 688**

**General**

**Safety:** This product is non-hazardous.  
**Storage:** 2-8°C. Do not freeze.  
**Catalogue Number:** ESC100  
**Preparation Date:** 5 December 2019  
**Expiration Date:** 3 April 2020  
**Volume:** 0.595 ml +/- 11 ul  
**Suspension media:** Buffered saline solution.  
**Sterilisation method:** Gamma irradiation.

| <b>Counts</b><br><small>(Method Ref: ESCG-001)</small> |            |            |
|--|------------|------------|
|  | Mean       | St.Dev.    |
| <b>Cryptosporidium count:</b>                          | <b>99</b>  | <b>2.2</b> |
| <b>DAPI staining:</b>                                  |            |            |
| <b>Cryptosporidium % +ve</b>                           | <b>100</b> | <b>%</b>   |

**Stock specifics**

**Organism:** *Cryptosporidium parvum*  
**Strain:** Iowa  
**Source:** Bovine  
**Shed date:** 20 November 2019  
**Purification method:** Discontinuous sucrose and cesium chloride centrifugation

**Quality Manager**



Lucy Millican  
 16 December 2019

**Water Industry Enquiries:**

Biopoint Pty Ltd  
 Unit 1 35-41 Waterloo Rd  
 Macquarie Park NSW 2113  
 Tel. 61 (0) 2 8877 9127  
**www.biopoint.com.au**  
 info@biopoint.com.au

