

**Cryptosporidium parvum**  
**10E8**

Product Information Sheet

**Batch number: B 804 - 506**

**Production Data**

**Volume:** 1ml  
**Viability:** Viable  
**Catalogue Number:** C10E8-V  
**Production Date:** 17 March 2020

**Stock specifics**

**Strain:** Iowa  
**Source:** Bovine  
**Date purified:** 20 January 2020  
**Date Excreted:** 20 January 2020  
**Purification Method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

**Quality Manager**



Lucy Millican  
27 March 2020

**Water Industry Enquiries:**

Biopoint Pty Ltd  
Unit 1 35-41 Waterloo Rd  
Macquarie Park NSW 2113  
Tel. 61 (0) 2 8877 9127  
**www.biopoint.com.au**  
info@biopoint.com.au

