

Giardia lamblia
10E6

Product Information Sheet

Batch number: B 813 - 510

Production Data

Volume: 1
Viability: Inactivated - Gamma Irradiation
Catalogue Number: G10E6
Production Date: 23 March 2020

Stock specifics

Strain: H3
Source: Gerbil
Date purified: 2 March 2020
Date Excreted: 2 March 2020
Purification Method: Sucrose and Percoll density gradient centrifugation

Quality Manager



Lucy Millican
27 March 2020

Water Industry Enquiries:

Biopoint Pty Ltd
Unit 1 35-41 Waterloo Rd
Macquarie Park NSW 2113
Tel. 61 (0) 2 8877 9127
www.biopoint.com.au
info@biopoint.com.au

