

Cryptosporidium parvum
10E7

Product Information Sheet

Batch number: B 814 - 509

Production Data

Volume: 1ml
Viability: Inactivated - Gamma Irradiation
Catalogue Number: C10E7
Production Date: 23 March 2020

Stock specifics

Strain: Iowa
Source: Bovine
Date purified: 4 March 2020
Date Excreted: 4 March 2020
Purification Method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Manager



Lucy Millican
27 March 2020

Water Industry Enquiries:

Biopoint Pty Ltd
Unit 1 35-41 Waterloo Rd
Macquarie Park NSW 2113
Tel. 61 (0) 2 8877 9127
www.biopoint.com.au
info@biopoint.com.au

