

Giardia lamblia
10E7

Product Information Sheet

Batch number: B 826 - 522

Production Data

Volume: 1.96ml
Viability: Inactivated - Gamma Irradiation
Catalogue Number: G10E7
Production Date: 21 July 2020

Stock specifics

Strain: H3
Source: Gerbil
Date purified: 3 July 2020
Date Excreted: 3 July 2020
Purification Method: Sucrose and Percoll density gradient centrifugation

Quality Manager



.....
Lucy Millican
30 July 2020

Manufactured By:
Biopoint Pty Ltd
Suite 16, 13A Narabang Way
Belrose, Sydney, NSW 2285
Tel: +61 2 8316 7939
www.biopoint.com.au
info@biopoint.com.au

