



Batch Number: **B 155**

**Fluorescence analysis - Mean Fluorescence Intensity**

<i>Cryptosporidium</i> oocysts:	<b>1269.41</b>
<i>Giardia</i> cysts:	<b>1007.73</b>

**Non-specific binding assay**

Number of particles detected:	<b>0</b>
Control number of particles detected:	<b>0</b>

**General**

Safety: This product contains sodium azide. Handle with care.

Storage: 2-8°C. Do not freeze.

Catalogue Number: EST-CG80-D, EST-CG20-D, JPN-EST-CG80-D, JPN-EST-CG20-D

Preparation Date: 3 May 2022

Expiration Date\*: 4 May 2023

*\*Expiration Date applies to all kit contents*

**Component Lots**

Mounting Medium Lot	MM-	1002
Wash Buffer Lot	WB-	1002
Postive Control Lot	PC-	1006
DAPI Lot	D-	1002

**Quality Manager**

Lucy Millican  
10 May 2022

**Manufactured By:**

Biopoint Pty Ltd  
Suite 16, 13A Narabang Way, Belrose, NSW, 2085  
Tel. +61 2 8316 7939  
[www.biopoint.com.au](http://www.biopoint.com.au)  
[info@biopoint.com.au](mailto:info@biopoint.com.au)