

Giardia lamblia
10E7

Product Information Sheet

Batch number: B 932 - 607

Production Data

Volume: 1.980ml
Viability: Inactivated - Gamma Irradiation
Catalogue Number: G10E7
Production Date: 1 August 2022

Stock specifics

Strain: H3
Source: Gerbil
Date purified: 15 July 2022
Date Excreted: 15 July 2022
Purification Method: Sucrose and Percoll density gradient centrifugation

Quality Manager



Lucy Millican
11 August 2022

Manufactured By:

Biopoint Pty Ltd
Suite 16, 13A Narabang Way
Belrose, Sydney, NSW 2085
Tel: +61 2 8316 7939
www.biopoint.com.au
info@biopoint.com.au



Issue:2 12 May 2020