

Cryptosporidium parvum
10E7

Product Information Sheet

Batch number: B 954 - 621

Production Data

Volume: 1ml
Viability: Viable
Catalogue Number: C10E7-V
Production Date: 21 November 2022

Stock specifics

Strain: Iowa
Source: Bovine
Date purified: 14 October 2022
Date Excreted: 11 October 2022
Purification Method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Manager



Lucy Millican
21 November 2022

Manufactured By:

Biopoint Pty Ltd
Suite 16, 13A Narabang Way
Belrose, Sydney, NSW 2085
Tel: +61 2 8316 7939
www.biopoint.com.au
info@biopoint.com.au



Issue:2 12 May 2020