

Giardia lamblia
10E6

Product Information Sheet

Batch number: B 964 - 626

Production Data

Volume: 1ml
Viability: Inactivated - Gamma Irradiation
Catalogue Number: G10E6
Production Date: 13 February 2023

Stock specifics

Strain: H3
Source: Gerbil
Date purified: 2 February 2023
Date Excreted: 1 February 2023
Purification Method: Sucrose and Percoll density gradient centrifugation

Quality Manager



Lucy Millican
21 February 2023

Manufactured By:

Biopoint Pty Ltd
Suite 16, 13A Narabang Way
Belrose, Sydney, NSW 2085
Tel: +61 2 8316 7939
www.biopoint.com.au
info@biopoint.com.au



Issue:2 12 May 2020