



Batch Number: **B 161**

Fluorescence analysis - Mean Fluorescence Intensity

| | |
|---------------------------------|----------------|
| <i>Cryptosporidium</i> oocysts: | 1198.70 |
| <i>Giardia</i> cysts: | 1192.43 |

Non-specific binding assay

| | |
|---------------------------------------|----------|
| Number of particles detected: | 0 |
| Control number of particles detected: | 0 |

General

Safety: This product contains sodium azide. Handle with care.

Storage: 2-8°C. Do not freeze.

Catalogue Number: EST-CG80-D, EST-CG20-D, JPN-EST-CG80-D, JPN-EST-CG20-D

Preparation Date: 6 June 2023

Expiration Date*: 6 June 2024

**Expiration Date applies to all kit contents*

Component Lots

| | | |
|---------------------|-----|------|
| Mounting Medium Lot | MM- | 1005 |
| Wash Buffer Lot | WB- | 1005 |
| Postive Control Lot | PC- | 1010 |
| DAPI Lot | D- | 1005 |

Quality Manager

Lucy Millican
9 June 2023

Manufactured By:

Biopoint Pty Ltd
Suite 16, 13A Narabang Way, Belrose, NSW, 2085
Tel. +61 2 8316 7939
www.biopoint.com.au
info@biopoint.com.au