

# **Cryptosporidium parvum** 10E7

## Product Information Sheet

**Batch number: B 991 - 639**

### Production Data

**Volume:** 1ml  
**Viability:** Inactivated - Gamma Irradiation  
**Catalogue Number:** C10E7  
**Production Date:** 31 July 2023

### Stock specifics

**Strain:** Iowa  
**Source:** Bovine  
**Date purified:** 26 June 2023  
**Date Excreted:** 26 June 2023  
**Purification Method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

### Quality Manager



Lucy Millican  
10 August 2023

### Manufactured By:

Biopoint Pty Ltd  
Suite 16, 13A Narabang Way  
Belrose, Sydney, NSW 2085  
Tel: +61 2 8316 7939  
[www.biopoint.com.au](http://www.biopoint.com.au)  
[info@biopoint.com.au](mailto:info@biopoint.com.au)



Issue:2 12 May 2020