

# **Cryptosporidium parvum** 10E7

## Product Information Sheet

**Batch number: B 993 - 642**

### **Production Data**

**Volume:** 1ml  
**Viability:** Inactivated - Gamma Irradiation  
**Catalogue Number:** C10E7  
**Production Date:** 11 September 2023

### **Stock specifics**

**Strain:** Iowa  
**Source:** Bovine  
**Date purified:** 26 July 2023  
**Date Excreted:** 26 July 2023  
**Purification Method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

### **Quality Manager**



Lucy Millican  
27 September 2023

### **Manufactured By:**

Biopoint Pty Ltd  
Suite 16, 13A Narabang Way  
Belrose, Sydney, NSW 2085  
Tel: +61 2 8316 7939  
**www.biopoint.com.au**  
info@biopoint.com.au



Issue:2 12 May 2020