

Cryptosporidium parvum
10E6

Product Information Sheet

Batch number: **B1001 - 645**

Production Data

Volume: 1ml
Viability: Viable
Catalogue Number: C10E6-V
Production Date: 23 October 2023

Stock specifics

Strain: Iowa
Source: Bovine
Date purified: 12 October 2023
Date Excreted: 12 October 2023
Purification Method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Manager



Lucy Millican
24 October 2023

Manufactured By:

Biopoint Pty Ltd
Suite 16, 13A Narabang Way
Belrose, Sydney, NSW 2085
Tel: +61 2 8316 7939
www.biopoint.com.au
info@biopoint.com.au



Issue:2 12 May 2020