

**Giardia lamblia**  
10E7

## Product Information Sheet

**Batch number: B1011 - 657**

### Production Data

**Volume:** 2.186ml  
**Viability:** Inactivated - Gamma Irradiation  
**Catalogue Number:** G10E7  
**Production Date:** 5 January 2024

### Stock specifics

**Strain:** H3  
**Source:** Gerbil  
**Date purified:** 1 December 2023  
**Date Excreted:** 1 December 2023  
**Purification Method:** Sucrose and Percoll density gradient centrifugation

### Quality Manager



Lucy Millican  
10 January 2024

### Manufactured By:

Biopoint Pty Ltd  
Suite 16, 13A Narabang Way  
Belrose, Sydney, NSW 2085  
Tel: +61 2 8316 7939  
**www.biopoint.com.au**  
info@biopoint.com.au



Issue:2 12 May 2020