

# **Cryptosporidium parvum** **C10E7**

## Product Information Sheet

**Batch number: B1036 - 671**

### **Production Data**

**Volume:** 1ml  
**Viability:** Inactivated - Gamma Irradiation  
**Catalogue Number:** C10E7  
**Production Date:** 17 June 2024

### **Stock specifics**

**Strain:** Iowa  
**Source:** Bovine  
**Date purified:** 29 May 2024  
**Date Excreted:** 29 May 2024  
**Purification Method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

### **Quality Manager**



Lucy Millican

25 June 2024

### **Manufactured By:**

Biopoint Pty Ltd  
Suite 16, 13A Narabang Way  
Belrose, Sydney, NSW 2085  
Tel: +61 2 8316 7939  
**[www.biopoint.com.au](http://www.biopoint.com.au)**  
[info@biopoint.com.au](mailto:info@biopoint.com.au)



Issue:2 12 May 2020