

Cryptosporidium parvum
C10E7

Product Information Sheet

Batch number: B1037 - 673

Production Data

Volume: 1ml
Viability: Inactivated - Gamma Irradiation
Catalogue Number: C10E7
Production Date: 29 July 2024

Stock specifics

Strain: Iowa
Source: Bovine
Date purified: 11 July 2024
Date Excreted: 11 July 2024
Purification Method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Manager



Lucy Millican

7 August 2024

Manufactured By:

Biopoint Pty Ltd
Suite 16, 13A Narabang Way
Belrose, Sydney, NSW 2085
Tel: +61 2 8316 7939
www.biopoint.com.au
info@biopoint.com.au



Issue:2 12 May 2020