

Giardia lamblia
G10E6

Product Information Sheet

Batch number: B1066 - 686

Production Data

Volume: 1ml
Viability: Inactivated - Gamma Irradiation
Catalogue Number: G10E6
Production Date: 18 February 2025

Stock specifics

Strain: H3
Source: Gerbil
Date purified: 6 February 2025
Date Excreted: 3 February 2025
Purification Method: Sucrose and Percoll density gradient centrifugation

Quality Manager



Lucy Millican
26 February 2025

Manufactured By:

Biopoint Pty Ltd
Suite 16, 13A Narabang Way
Belrose, Sydney, NSW 2085
Tel: +61 2 8316 7939
www.biopoint.com.au
info@biopoint.com.au



Issue:2 12 May 2020