

Giardia lamblia
G10E7

Product Information Sheet

Batch number: B1070 - 689

Production Data

Volume: 3.4mL
Viability: Inactivated - Heat Inactivation
Catalogue Number: G10E7
Production Date: 24 February 2025

Stock specifics

Strain: H3
Source: Gerbil
Date purified: 1 February 2024
Date Excreted: 31 January 2024
Purification Method: Sucrose and Percoll density gradient centrifugation

Quality Manager



Lucy Millican
27 February 2025

Manufactured By:

Biopoint Pty Ltd
Suite 16, 13A Narabang Way
Belrose, Sydney, NSW 2085
Tel: +61 2 8316 7939
www.biopoint.com.au
info@biopoint.com.au



Issue:2 12 May 2020